



positioning
patients for life.®

MOLDCARE® Cushions

Provides Customized Patient Support

MOLDCARE® Cushions provide customized head and neck support. Patient comfort is a critical component of patient positioning. The MOLDCARE® cushion is composed of a soft fabric that contains expanded polystyrene beads coated in a moisture-cured resin. When activated with water, the cushion becomes rigid, conforming to the contours of the patient's head and neck. It can be used alone or on a Silverman Head Support to create the most comfortable and repeatable head and neck patient specific support.



RT-4492S
Small, 15 cm x 20 cm



RT-4492
20 cm x 25 cm



RT-4492U
20 cm x 35 cm
(Shown molded & with Q2 Headrest)



RT-4492T
Head & Shoulder,
Trapezoid, 58 cm x 65 cm



RT-4492C
Head & Shoulder S-Frame

RT-4492-B
Spray Bottle



RT-4492	MOLDCARE® Head Cushion, 20 cm x 25 cm
RT-4492M	MOLDCARE® Head Cushion, 20 cm x 45 cm
RT-4492S	MOLDCARE® Head Cushion, Small, 15 cm x 20 cm
RT-4492T	MOLDCARE® Trapezoid Head & Shoulder Cushion, 58 cm x 65 cm
RT-4492L	MOLDCARE® Cushion, Large, 40 cm x 60 cm
RT-4492XL	MOLDCARE® Cushion, Extra-Large, 60 cm x 60 cm
RT-4492U	MOLDCARE® U Head Cushion, 20 cm x 35 cm
RT-4492C	MOLDCARE® Head and Shoulder S-Frame Cushion
RT-4492-B	Spray Bottle for MOLDCARE®



Silverman Head Support

Provides customized patient support

Silverman head and neck supports are available in a variety of sizes to fit your positioning needs.

- The Silverman Q1 and Q2 minimize overstretching of the thermoplastic mask while simultaneously offering cranial support
- Silverman & MOLDCARE® together are lower attenuating than most standard foam headrests
- Q Headrests feature a posterior cranial indexing notch for enhanced reproducibility of MOLDCARE® cushions
- Q Headrests extend to the top of the cranium for added support

Silverman Q Head Supports



RT-4491-Q1
Silverman Q1



RT-4491-Q2
Silverman Q2



RT-4491R-M*
Silverman 'M' Low-Profile
Head Support



RT-4491-Q3-2G
Silverman Q3



RT-4491-Q4-2G
Silverman Q4



*Exclusive low-profile design allows the occiput to rest below the top plane of the board. By lowering the patient's head, the mask is stretched a shorter distance during the molding process. Less stretch means a more comfortable mask and an overall better fit for your patient.

Silverman Head Support

RT-4491R2 Silverman Head Supports, set of 6, sizes A-F



RT-4491R2-A
size A



RT-4491R2-B
size B



RT-4491R2-C
size C



RT-4491R2-D
size D



RT-4491R2-E
size E



RT-4491R2-F
size F