



positioning  
patients for life.®

# Stradivarius™

for SBRT

***Complete Solution for  
Hypofractionated Non-Coplanar Treatment***




# Stradivarius™



## Variable Height Bridges

### RT-4558KV2 & RT-4558KV3

- Bridge Dimensions\*:
  - 38 cm (RT-4558KV2) min: 28cm, max: 38cm
  - 48 cm (RT-4558KV3) min: 38cm, max: 48cm
- \*Measurements are from underside of bridge to the top of the kVue*
- Adjustable Height
- Easy Single Motion Lever
- Horizontal & Vertical Laser Alignment Scales
- Low Attenuation to Enable Treatment Through Bridge




## Compression Paddle

### RT-4558KV4 (shown)

### RT-4558KV4-01

*Low profile version for use in a small bore CT*

- Safety Release Button
- Tactile Compression
- Incremented Lateral Placement Along Bridge







## Thigh Compression VacQfix™ Cushion

RT-4517-4079VF15

- Complete Leg Immobilization
- Adjustable Along Bridge
- One-Step™ Valve
- Snap on Attachment



## Whole Body VacQfix™ Cushion

RT-4517-230120T80F90


- Complete Body Immobilization
- Flared Top for Arm Immobilization
- One-Step™ Valve
- Range of Sizes



## kVue™ Stradivarius SBRT Couch Top Insert

RT-4558KV1

- Compatible with kVue™ IGRT Couch Top
- Indexing Every 70 mm- Double Your Indexing Options
- Laser Alignment Scales
- Central Cut-Out for Unobstructed Treatment Paths including Non-Coplanar Fields



Patented

## Stradivarius™ SBRT Overlay

RT-4558-1

- Overlays any couch top that uses Two-Pin Indexing Bars
- Indexing Every 70 mm - Double Your Indexing Options
- Laser Alignment Scales
- Central Cut-Out for Unobstructed Treatment Paths including Non-Coplanar Fields



# Stradivarius™



## ZiFix™, Abdominal Thoracic Motion Control System

### RT-4558CB02

- Suppresses motion of lungs and liver for imaging & radiotherapy
- Compact, one-body design
- Three belt sizes



## SofTouch™ Knee Wedges (large & small) & Foot Block

### RT-4489KW (large)

### RT-4489KW-01 (small)

### & RT-4489FB

- Ergonomic design
- Indexes to locating bar in 35 mm increments
- Clip-on leg straps for added immobilization (included)
- Centering groove to aid in laser alignment



## ArmShuttle™ Elite

### RT-4546C


- Lightweight
- Ergonomic design
- MR Safe & CT compatible



## VacQfix™ Indexer

### RT-4517-IND01

- Quick, easy vacuum cushion indexing
- Adaptable for a variety of setups



## Leg Roll Positioner Cushions

### RT-4474-6 (60 cm x 15 cm, 24 in x 6 in)

### RT-4474-8 (60 cm x 20 cm, 24 in x 8 in)

- Durable for everyday use
- Maintain consistent shape for repeatability
- Help flatten the lower lumbar region